Pierre-François Duc & Benjamin Schmidt

LabGui

or automating measurements with Python



Manual measurements



LabView



Python

Automating your experiment with LabGui

May 19, Tuesday, noon


3600 University, Rutherford Physics Building, room 103

We have funds to help you cover transportation costs. If interested contact us at:

rtech@physics.mcgill.ca


https://www.facebook.com/RTechMcGill

Determine resistance




Start measurement #1






Set current value I to 1A







Measure voltage V






Plot your first point







Start measurement #2






Set current value I to 2A







Measure voltage V

#		V
1	1.0	10.3
2	2.0	21.3



Plot your second point






How many measurements you say one needs to do to reduce statistical uncertainty?


so, continue...


And here is your data

#		V
1	1.0	10.3
2	2.0	21.3
3	3.0	29.9
4	4.0	39.5
5	5.0	51.2



And here is your plot






Linear regression V = IR

$$y = X\beta$$


$$\mathbf{y} = \begin{pmatrix} y_1 \\ y_2 \\ \vdots \\ y_n \end{pmatrix}, \quad \mathbf{X} = \begin{pmatrix} \mathbf{x}_1^{\mathrm{T}} \\ \mathbf{x}_2^{\mathrm{T}} \\ \vdots \\ \mathbf{x}_n^{\mathrm{T}} \end{pmatrix} = \begin{pmatrix} x_{11} & \cdots & x_{1p} \\ x_{21} & \cdots & x_{2p} \\ \vdots & \ddots & \vdots \\ x_{n1} & \cdots & x_{np} \end{pmatrix}, \quad \boldsymbol{\beta} = \begin{pmatrix} \beta_1 \\ \beta_2 \\ \vdots \\ \beta_p \end{pmatrix}$$








And... Here you are!



V = I R

R = 10.12

Is it a job for a student?



Only two options left







Python / LabGui

Communicate with the device


Ethernet

GPIB









RS232



USB




Ohm's Law example

LabView or LabGui

Let's see how LabGui performs

Set & Measure with a Click




```
#empty array which will contain the measurement results
measured_voltages=[]

#create instrument objects using the communication port
current_source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")

#This is the experiment
I=1.0
current_source.set_current(I)
V=voltmeter.measure_voltage()
measured_voltages.append(V)
```

#code to perform one measurement

Set & Measure with a Click



Is that hard to automate?

Is that hard to automate?


Nope, with **Python** it's easy!

Just add a FOR loop.

```
#empty array which will contain the measurement results
measured_voltages=[]
#Value of the current in amperes
currents=[1,2,3,4,5]
#create instrument objects using the communication port
current_source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")
#This is the experiment
for I in currents:
    current_source.set_current(I)
    V=voltmeter.measure_voltage()
    measured_voltages.append(V)
    #plots the graph of V versus I
    plot(I,measured_voltages)
```

#code to perform five measurements


5 experimental points at once




What about fitting?

```
#empty array which will contain the measurement results
measured voltages=[]
#Value of the current in amperes
currents=[1,2,3,4,5]
#create instrument objects using the communication port
current_source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")
def linear(x,a,b):
                                                            #code added
    return a*x+b
#This is the experiment
for I in currents:
    current source.set current(I)
    V=voltmeter.measure_voltage()
    measured_voltages.append(V)
    #plots the graph of V versus I
    plot(I,measured voltages)
    #fit the function
    if fit button on==True:
        perform_fit(I,measured_voltages,name_of_function=linear)
                                                            #code added
```

Calculate and plot fit



Find fitting parameters




Ok, what about non-linear fitting?

Ok, what about **non-linear** fitting? Any function!


This is **Python**YOU control everything

```
#empty array which will contain the measurement results
measured_voltages=[]
#Value of the current in amperes
currents=[1,2,3,4,5]
#create instrument objects using the communication port
current_source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")
def non linear whatever(x,a,b):
                                            #code changed here
    return a*exp(b*x)
#This is the experiment
for I in currents:
    current source.set current(I)
    V=voltmeter.measure_voltage()
    measured_voltages.append(V)
    #plots the graph of V versus I
    plot(I,measured_voltages)
   #fit the function
    if fit button on==True:
       perform_fit(I,measured_voltages,name_of_function=non_linear_whatever)
                                            #code changed here
```


Choose a fitting function




Choose a fitting function




Find fitting parameters



Choose another fitting function



Calculate and plot new fit



Want more statistics?

Want more statistics?

Just add another **FOR** loop

```
repetitions=[1,2,3,4,5,6,7,8,9,10]
#empty array which will contain the measurement results
measured voltages=[]
#Value of the current in amperes
currents=[1,2,3,4,5]
#create instrument objects using the communication port
current_source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")
#This is the experiment
for j in repetitions:
                                                #code changed on
    for I in currents:
        current source.set current(I)
        V=voltmeter.measure voltage()
        measured_voltages.append(V)
        #plots the graph of V versus I
        plot(I, measured voltages)
```

Want more statistics?


Just add another **FOR** loop

Bonus:


Play with fits of already taken data while

experiment goes on in the bckg


Meanwhile in the background



Meanwhile in the background



Go grab a coffee while your experiment is going on!



Go grab a coffee while your experiment is going on!

Another bonus: LabGui will send you an email, if something happens!



repetitions=[1,2,3,4,5,6,7,8,9,10]


```
#empty array which will contain the measurement results
measured voltages=[]
#Value of the current in amperes
currents=[1,2,3,4,5]
#create instrument objects using the communication port
current source=CurrentSource("COM1")
voltmeter=Voltmeter("COM2")
#This is the experiment
for j in repetitions:
    for I in currents:
       current source.set current(I)
        V=voltmeter.measure voltage()
       measured voltages.append(V)
        #plots the graph of V versus I
        plot(I, measured voltages)
#Call a function which send you an email
send_email("you@some_server.org", message="Come back, the measurement is over ;)")
                           #all the code you need to send an email
```

Ohm's Law example

Why should I use LabGui?


Let's compare the key features

Compare LabView & LabGui



#1 ALL devices speak ONE language

Compare LabView & LabGui



#2 Project's complexity

#2 Project's complexity


Source of complications: PL vs GPL

LabGui = Programming Language LabView = Graphical Programming Language

LabGui

Code →

```
85 class LabGui(QtGui.QMainWindow):
86
       #The command window
87
       cmdwin=None
88
89
       outputfile=None
       def __init__(self):
90
91
           # run the initializer of the class inherited from6
92
           super(LabGui, self).__init__()
93
94
           self.settings = QSettings(self)
95
           self.settings.setValue("state", self.saveState())
96
97
           #defines the zone in which you can create widgets
98
           self.zoneCentrale = OtGui.OMdiArea()
```




Log X Log Y Log Y

Interface


LabView

Code →

Interface



LabView



Actually that diagram looks benign compared to some of the stuff that we use in our lab.

"LabView makes:

- the easy things easier and
- the hard things harder."

"LabView makes:

- the easy things easier and
- the hard things harder."
- The learning curve of GPL is lower, but so is the productivity

"LabView makes:


- the easy things easier and
- the hard things harder."
- The learning curve of GPL is lower, but so is the productivity
- GPL has less convenient and standard VERSION CONTROL

"LabView makes:

- the easy things easier and
- the hard things harder."
- The learning curve of GPL is lower, but so is the productivity
- GPL has less convenient and standard VERSION CONTROL

 It's cumbersome and slow to use GPL over SSH and/or on old computers

Compare LabView & LabGui



#3

Proprietary vs Open Source software

 Full control of the program: you know, what is inside






LabGui

- Full control of the program: you know, what is inside
- It's free






LabGui





Use any text editor to modify your code






Use any text editor to modify your code no need to fight with your colleagues for the computer, where the latest version of LabView is installed

LabView compatibility issues

Table 1: Possible save/open combinations

Save in Version:	Open in Version:														
	5.0.x	5.1.x	6.0.x	6.1	7.0	7.1.x	8.0.x	8.2.x	8.5.x	8.6.x	2009	2010	2011	2012	2013
2.x	C	С	C	С	C	С	С	С	C+I	C+I	C+I	C+I	C+I	C+I	C+I
3.x	C	C	C	C	C	С	C	C	C+I	C+I	C+I	C+I	C+I	C+I	C+I
4.x	1	1	1	1	1	1	1	1	1	-1	1	1	1	- 1	1
5.0.x	1	1	1	1	1	1	1	1	1	1	1	1	1	1	1
5.1.x	S	1	1	1	1	1	1	1	L	1	1	1	1	- 1	1
6.0.x	M	S	1	1	1	1	1	1	1	1	1	1	1	1	1
6.1	M	M	S	1	1	1	1	1	1	1	1	1	1	1	1
7.0	М	M	М	S	1	1	1	1	1	1	1	1	1	1	1
7.1.x	M	M	M	M	S	1	1	1	1	1	1	1	1	1	1
8.0.x	M	М	M	M	M	S	1	1	1	1	1	1	1	1	1
8.2.x	M	M	M	M	M	M	S	1	1	1	1	1	1	1	1
8.5.x	M	M	М	M	M	M	S	S	1	1	1	1	1	1	1
8.6.x	M	M	M	M	M	M	S	S	S	1	1	1	1	1	1
2009	M	М	M	M	M	M	S	S	S	S	1	1	1	1	1
2010	M	M	M	M	M	M	S	S	S	S	S	1	1	1	1
2011	M	M	М	M	М	M	S	S	S	S	S	S	1	1	1
2012	M	M	M	M	M	M	S	S	S	S	S	S	S	1	1
2013	M	М	М	M	M	M	S	S	S	S	S	S	S	S	1


LabView compatibility issues



2013

Use all the power and support:

HUGE Python community





Use all the power and support:


- HUGE Python community
- Growing scientific community

Use all the power and support:

- HUGE Python community
- Growing scientific community


List of papers that used LabGui

- B. A. Schmidt, K. Bennaceur, S. Bilodeau, K. W. West, L. N. Pfeiffer, G. Gervais, "5/2 Fractional Quantum Hall Effect in the Corbino Geometry", http://arxiv.org/abs/1503.07775
- P-F Duc, M.Savard, M. Petrescu, B. Rosenow, A. Del Maestro, and G. Gervais, "Critical Flow and Dissipation in a Quasi-One-Dimensional Superfluid", http://arxiv.org/abs/1412.5124, accepted in Science Advances
- V. Tayari, N. Hemsworth, I. Fakih, A. Favron, E. Gaufres, G. Gervais, R. Martel, T. Szkopek, "Two-Dimensional Magnetotransport in a Black Phosphorus Naked Quantum Well", http://arxiv.org/abs/ 1412.0259, accepted in Nature Communications



Summary LabView & LabGui

Summary LabView & LabGui



Download LabGui

https://bitbucket.org/ RTechMcGill/ labgui/ src

Automating your experiment with LabGui

May 19, Tuesday, noon

3600 University, Rutherford Physics Building, room 103

We have funds to help you cover transportation costs. If interested contact us at:

rtech@physics.mcgill.ca

https://www.facebook.com/RTechMcGill

www.hep.physics.mcgill.ca/RTech facebook.com/RTechMcGill



Video created by: Pierre-François Duc & Igor Kozlov

> with contributions by: Julien Lhermitte Anna Mkrtchyan Hélène Seiler